

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Biphenylacetic acid - 98%

Product Number: 196487
CAS Number: 5728-52-9
MDL: MFCD00004351
Formula: C₁₄H₁₂O₂
Formula Weight: 212.24 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals
Appearance (Form)	Powder or Crystals
Purity (Titration by NaOH)	97.5 - 102.5 %

Specification: PRD.0.ZQ5.10000043510

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.