


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Triethyloxonium tetrafluoroborate solution - 1.0 M in methylene chloride

Product Number: 176230



MDL: MFCD00044423

Formula: C6H15BF4O

Formula Weight: 189.99 g/mol

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Very Dark Yellow and Very Faint Brown-Yellow to Brown-Yellow to Brown
Appearance (Form)	Liquid
Proton NMR spectra	Conforms to Structure
Purity Titration by NAOH	0.95 - 1.10 M
Expiration Date Period	-----
9 months	

Specification: PRD.6.ZQ5.10000011949

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

