

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Sheep IgG (whole molecule)-Peroxidase antibody produced in donkey - affinity isolated antibody, buffered aqueous solution

Product Number: A3415

MDL: MFCD00162847

Storage Temperature: -20 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Light Brown
Appearance (Form)	Liquid
Protein Content	
Conjugate concentration prior to the addition of BSA	
Report Result	
Direct ELISA Titer-Sheep IgG	10000 - 30000
1:10,000 - 1:30,000 working dilution versus 5 µg/ml sheep IgG.	

Specification: PRD.4.ZQ5.10000006152

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

