


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methanol-d₃ - 99.8 atom % D

Product Number: 343854
CAS Number: 1849-29-2
MDL: MFCD00055877
Formula: CD₃HO
Formula Weight: 35.06 g/mol



| TEST | Specification |
|------------------------|---------------|
| Appearance (Clarity) | Clear |
| Appearance (Color) | Colorless |
| Appearance (Form) | Liquid |
| Purity (GC) | ≥ 99 % |
| Proton NMR Spectrum | Conforms |
| Deuterium NMR Spectrum | Conforms |
| Proton NMR | ≥ 99.8 % |
| Water Content | ≤ 500 ppm |

Specification: PRD.0.ZQ5.10000065213

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.