

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cadmium tungstate - -325 mesh

Product Number: 401129
CAS Number: 7790-85-4
MDL: MFCD00016002
Formula: CdO4W
Formula Weight: 360.25 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Faint Green to Green to Brown
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
ICP Master Analysis	Confirmed
Confirms Cadmium and Tungsten Components	
Loss on Ignition	
1100 Degrees Celsius	
Particle Size: -325 Mesh	Conforms

Specification: PRD.1.ZQ5.1000067197

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.