


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Trichloroacetic acid solution - 6.1 N

Product Number: T0699

Formula: C₂HCl₃O₂
Formula Weight: 163.39 g/mol

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Normality by Titration by Sodium Hydroxide	5.8 - 6.4 Normal
Concentration	0.95 - 1.05 g/ml

Specification: PRD.0.ZQ5.10000092498

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.