


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2-Butanone - for HPLC, ≥99.7%

Product Number: 34861  
CAS Number: 78-93-3  
MDL: MFCD00011648  
Formula: C<sub>4</sub>H<sub>8</sub>O  
Formula Weight: 72.11 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
UV Absorbance 350nm	≤ 0.01
UV Absorbance 340nm	≤ 0.10
UV Absorbance 335nm	≤ 0.30
Purity (GC)	≥ 99.70 %
Residue on Evaporation	≤ 0.001 %
Water (by Karl Fischer)	≤ 0.1 %
Free Acid (as Acetic Acid)	≤ 0.001 %
as Acetic Acid	
Recommended Retest Period	-----
42 Months	

Specification: PRD.3.ZQ5.10000022523

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

