

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

3,3',5-Triiodo-L-thyronine sodium salt -  $\gamma$ -irradiated, powder, BioXtra, suitable for cell culture

Product Number:

**T5516**

CAS Number:

55-06-1

MDL:

MFCD00002594

Formula:


C<sub>15</sub>H<sub>11</sub>I<sub>3</sub>NNaO<sub>4</sub>

Formula Weight:

672.96 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

White to Off White and Pink

Appearance (Form)

Powder

Solubility (Turbidity)

Clear

1H NMR Spectrum

Conforms to Structure

Performance Testing

Pass

Sterilized (Gamma Irradiated)

Pass

(Irradiation dose 25.0 to 35.0 kGy)

Purity (HPLC)

≥ 90 %

Specification: PRD.1.ZQ5.10000062476

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.