

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Rabbit IgG (whole molecule)-FITC antibody produced in goat - IgG fraction of antiserum, buffered aqueous solution

Product Number: **F6005**

MDL: MFCD00162786

Storage Temperature: -20 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Faint Yellow to Yellow
Appearance (Form)	Solution
F/P Molar Ratio	3.0 - 5.0
A280/A496 UV Absorption	1.0 - 1.5
Protein Concentration by E280	10 - 20 mg/ml
Immunofluorescence (Indirect) using Mouse Spleen Cells	Title Minimum Dilution 1:40

Specification: PRD.1.ZQ5.10000010938

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.