


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Z-Phe-Arg 7-amido-4-methylcoumarin hydrochloride - kallikrein substrate

Product Number: C9521
CAS Number: 70382-26-2
MDL: MFCD00077030
Formula: C₃₃H₃₆N₆O₆ · HCl
Formula Weight: 649.14 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
% Water	≤ 10.0
Solubility (Color)	Colorless
Solubility (Turbidity) 20 mg/mL in MeOH	Clear
Carbon (anhydrous)	59.5 - 63.8 %
Nitrogen (anhydrous)	12.2 - 14.2 %
Purity (HPLC)	≥ 95 %

Specification: PRD.1.ZQ5.10000042712

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

