

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

ProteoMas™ MALDI Calibration Kit for Exactive Series MS, LTQ XL and LTQ Hybrids – MALDI MS Calibration Standard

Product Number: MSCAL4

Formula Weight: 1,046.2 g/mol

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Suitability	Suitable
Suitable for use in MALDI Mass Spectrometry. vMALDI instrument was successfully calibrated using the protocol and reagents supplied in the kit	
Note	Conforms
Exactive is a trademark of Thermo Fisher Scientific, Inc..	

Specification: PRD.1.ZQ5.10000005966

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.