


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Perfluorodecanoic acid - 98%

Product Number: 177741
CAS Number: 335-76-2
MDL: MFCD00004175
Formula: C₁₀HF₁₉O₂
Formula Weight: 514.08 g/mol



TEST

Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals or Chunk(s) or Powder with Chunk(s) or Crystals with Chunk(s)
Infrared spectrum	Conforms to Structure
Titration with NaOH	97.5 - 102.5 %
Purity (GC)	≥ 97.5 %
Solubility (Turbidity) 10% In Methanol	Clear to Slightly Hazy
Solubility (Color)	Colorless to Faint Yellow

Specification: PRD.1.ZQ5.10000084041

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

