

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Chaperonin 60 from Escherichia coli - &gt;95% (SDS-PAGE), recombinant, expressed in E. coli overproducing strain, lyophilized powder

Product Number:

**C7688**

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White to White w/ Yellow Cast
Gel Electroph. (Gel)	SDS-Page
Gel Electroph. (Stain)	Coomassie
Gel Electroph. (Purity)	≥ 95.00 %
Vial Content (Size)	1.000 - 1.000
Vial Content (Unit Of Measure)	mg
Vial Content (Actual Content)	1.000 - 1.300
Vial Content (Testing Method)	Bradford
Suitability (Result)	Pass
Suitability (Test)	Refolding (w/GroES)
Recommended Retest Period	-----
2 Years	

Specification: PRE.2.ZQ5.10000003807

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.