

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-FLIP α , C-Terminal antibody produced in rabbit - IgG fraction of antiserum, buffered aqueous solutionProduct Number: **F6550**

MDL: MFCD02096351

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Solution in PBS containing 0.02% Sodium Azide

Protein Content	Conforms
-----------------	----------

1 mg/mL

Working Dilution	Conforms
------------------	----------

by Direct Immunoblotting

0.5-1.0 ug/mL (or 1:500-1:2000) using total HELA cell lysate

Specification: PRD.0.ZQ5.10000076468

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.