


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
(Hydroxypropyl)methyl cellulose - average  $M_n \sim 90,000$

Product Number: **423181**  
CAS Number: 9004-65-3  
MDL: MFCD00131360



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Viscosity	13275 - 24780 cps
At 20 Degrees Celsius	

Specification: PRD.0.ZQ5.10000100072

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.