

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

2-Ethylhexyl acrylate - 98%, contains  $\geq 0.001$ - $\leq 0.11\%$  monomethyl ether hydroquinone as stabilizer

Product Number:

290815

CAS Number:

103-11-7

MDL:


MFCD00009495

Formula:

C<sub>11</sub>H<sub>20</sub>O<sub>2</sub>

Formula Weight:

184.28 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

 $\geq 97.5\%$ 

Stabilizer

0.001 - 0.110 %

Monomethyl Ether Hydroquinone (MEHQ)

Specification: PRD.0.ZQ5.10000012833

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.