

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Anisotropic Silicon Etchant

Product Number: 667471

TEST	Specification
Appearance (Color)	Conforms to Requirements
Colorless to Yellow to Red to Amber	
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Miscellaneous Supplier Data	Conforms
Ethylenediamine 57.5% - 58.5%	
Miscellaneous Supplier Data	Conforms
Pyrocatechol 20% - 21%	
Miscellaneous Supplier Data	Conforms
Pyrazine < 0.5%	
Water 20% - 21%	Conforms
Trace Metal Analysis	
Report Value	

Specification: PRD.0.ZQ5.10000061559

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.