

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N-[Tris(hydroxymethyl)methyl]acrylamide - contains $\leq 7\%$ KCl, 93%

Product Number: 364959
CAS Number: 13880-05-2
MDL: MFCD00081045
Formula: C₇H₁₃NO₄
Formula Weight: 175.18 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	43.2 - 51.8 %
Nitrogen	7.3 - 8.7 %

Specification: PRD.0.ZQ5.10000027952

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.