

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cadmium - granular, 30-80 mesh, ≥99%

Product Number: 414891
CAS Number: 7440-43-9
MDL: MFCD00010914
Formula: Cd
Formula Weight: 112.41 g/mol

Cd

TEST	Specification
Appearance (Color) Silver and/or Grey	Conforms to Requirements
Appearance (Form) Metallic Granules	Conforms to Requirements
Particle Size 30-80 Mesh	Confirmed
ICP Major Analysis Confirms Cadmium Component	Confirmed
Purity > = 99% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 10000.0 ppm

Specification: PRD.2.ZQ5.10000039715

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.