


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Azidomethyl phenyl sulfide - 95%

Product Number: 244546
CAS Number: 77422-70-9
MDL: MFCD00009700
Formula: C7H7N3S
Formula Weight: 165.22 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Brown
Appearance (Form)	Liquid
Powder or Crystals	
Infrared spectrum	Conforms to Structure
Carbon	47.8 - 53.9 %
Nitrogen	23.9 - 27.0 %
Sulfur	18.2 - 20.6 %
Purity (HPLC)	≥ 94.0 %

Specification: PRD.1.ZQ5.10000044648

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.