


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calcium nitride – –200 mesh, 95%

Product Number: 415103
CAS Number: 12013-82-0
MDL: MFCD00015985
Formula: Ca₃N₂
Formula Weight: 148.25 g/mol



TEST	Specification
Appearance (Color)	Brown
Appearance (Form)	Powder
Compleximetric EDTA % Ca	76.7 - 86.0 %
X-Ray Diffraction	Conforms to Structure
Size -200 mesh	Confirmed

Specification: PRD.1.ZQ5.10000025877

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.