

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Lauryl methacrylate - contains 500 ppm MEHQ as inhibitor, 96%

**Product Number:****291811**

CAS Number:

142-90-5

MDL:


MFCD00008972

Formula:

C<sub>16</sub>H<sub>30</sub>O<sub>2</sub>

Formula Weight:

254.41 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 95.5 %

MEHQ as Inhibitor

450 - 550 ppm

Specification: PRD.1.ZQ5.10000008736

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.