

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

D-Glucuronic acid - ≥98% (GC)

 Product Number:
 G5269

 CAS Number:
 6556-12-3

 MDL:
 MFCD00077778

Formula: C6H10O7 Formula Weight: 194.14 g/mol OH OH

TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)
100 mg/mL, H2O
Specific Rotation

c = 6, H2O, equilibrium value

Purity (GC)

1H NMR Spectrum

Recommended Retest Period

3 years

Specification: PRD.2.ZQ5.10000007972

Specification

White to Off White Powder Colorless to Faint Yellow Clear to Slightly Hazy

35 - 37 °

> 98 %

Conforms to Structure

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.