

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethyl 2-trans-4-cis-decadienoate - natural, ≥95%, FG

Product Number: W314811
CAS Number: 3025-30-7
MDL: MFCD00082225
Formula: C12H20O2
Formula Weight: 196.29 g/mol



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.482 - 1.488
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 95.0 %
Trans/cis isomer	
Organoleptics	Pass
Color Test	≤ 500 APHA
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.2.ZQ5.10000009909

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.