

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Pentaerythritol tetrakis(2-bromoisobutyrate) - 97%

Product Number: **723193**  
CAS Number: 243991-62-0

Formula: C<sub>21</sub>H<sub>32</sub>Br<sub>4</sub>O<sub>8</sub>  
Formula Weight: 732.09 g/mol



### TEST

### Specification

Appearance (Color)

Conforms to Requirements

White to Off-White to Yellow

Appearance (Form)

Powder to Crystalline Powder

Infrared spectrum

Conforms to Structure

Titration by AgNO<sub>3</sub>

96.5 - 103.5 %

After Oxygen Combustion

Purity (HPLC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000012887

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.