


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium arsenate monobasic

Product Number: A6631
CAS Number: 7784-41-0
MDL: MFCD00036296
Formula: H₂AsKO₄
Formula Weight: 180.03 g/mol



TEST	Specification
------	---------------

Appearance (Color) White to Off-White	Conforms
Appearance (Form)	Powder
% Arsenic	40.8 - 42.5

Specification: PRD.1.ZQ5.10000030792

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.