

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
α-Casein Dephosphorylated from bovine milk - lyophilized powder

Product Number: **C8032**

MDL: MFCD00130477

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear to Hazy
20 mg/ml of 0.1 M TRIS, pH 7.5	
Protein by Biuret	≥ 70 %
Phosphorus (P)	
Dephosphorylated	≥ 80 %

Specification: PRD.0.ZQ5.10000103693

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

