

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
17-Azido-3,6,9,12,15-pentaoxaheptadecanoic acid - ≥95%

Product Number: **910015**
CAS Number: 217180-81-9



Formula: C₁₂H₂₃N₃O₇
Formula Weight: 321.33 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid or Semi-Solid or Paste or Solid
¹ H NMR Spectrum	Conforms to Structure
¹³ C NMR Spectrum	Conforms to Structure
Purity (HPLC)	≥ 95.0 %
Titration with NaOH	95.0 - 105.0 %
LCMS	Confirmed
Confirms M-1	

Specification: PRD.0.ZQ5.10000101391

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.