

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Amphotericin B solution – 250 µg/mL in deionized water, 0.1 µm filtered, BioReagent, suitable for cell culture

Product Number:

A2942

Formula:

C₄₇H₇₃NO₁₇

Formula Weight:

924,08 g/mol

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Turbidity)	Clear
Appearance (Form)	Solution
Appearance (Color)	Yellow
pH	9.0 - 12.0
Potency	250 - 300 µg/ml
Endotoxin Level at 10 ml/L	≤ 0.5 EU/ml
Cell Culture Test NLT 50% of control	≥ 50 %
Cell Line	Attached Cell Line
Expiration Date Period 18 months	-----

Specification: PRE.0.ZQ5.10000020661

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.