

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)**Product Specification**

Product Name:

Potassium hexachloroiridate(III)

**Product Number:****342262**

CAS Number:

14024-41-0

MDL:

MFCD00038699

Formula:

Cl<sub>6</sub>IrK<sub>3</sub>

Formula Weight:

522.23 g/mol



<b>TEST</b>	<b>Specification</b>
Appearance (Color) Green to Brown	Conforms to Requirements
Appearance (Form) Powder or Crystals	Conforms to Requirements
Gravimetric Analysis % Ir After Reduction With Magnesium	35.5- 38.1 %
ICP Major Analysis Confirms K & Ir Components	Confirmed

**Remarks:**

Specification Date : 12/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.