


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Aminocarbonylphenylboronic acid

Product Number: 683868
CAS Number: 351422-73-6
Formula: C₇H₈BNO₃
Formula Weight: 164.95 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

White to Yellow and Faint Beige to Beige
and Faint Brown to Light Brown

Appearance (Form)

Crystal to Powder

Proton NMR spectrum

Conforms to Structure

Purity (Titration by NaOH)

≥ 95.0 %

Note

Note

May Contain Varying Amounts of Anhydride

Specification: PRD.0.ZQ5.10000038077

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.