

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Bromo-4-chloro-3-indolyl β-D-glucuronide sodium salt - tablet

Product Number: B8174
CAS Number: 129541-41-9
MDL: MFCD00135782
Formula: C₁₄H₁₂BrClNaO₇
Formula Weight: 444.59 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Round Tablet
3/16"	
Total Tablet Weight	36 - 44 mg
Dissolution Time	< 10 min
One Tablet Dissolved in 5 mL DH ₂ O	
Suitability	Pass
Suitable for detection of Gus Gene expression; compares to untableted substrate.	
One tablet dissolved in 0.5 mL DH ₂ O	
Recommended Retest Period	-----
4 Years	

Specification: PRD.0.ZQ5.10000031353

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.