

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Fluorophenol - 99%

Product Number: **F13207**
 CAS Number: 371-41-5
 MDL: MFCD00002316
 Formula: C₆H₅FO
 Formula Weight: 112.1 g/mol



TEST

Specification

Appearance (Color)	White to Yellow
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Solid, Crystalline Powder or Crystalline Granules and/or	Chunks
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
5% in H ₂ O	
Solubility (Color)	Conforms to Requirements

Colorless to Light Yellow or Light Yellow-Green

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.