

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
CARBON NANOTUBE SHEET, ALIGNED, SIZE 10&

Product Number: **901082**
CAS Number: 308068-56-6

Formula: C

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Sheets
Resistance	Confirmed
Sheet (ohm/sq) 10-10 ³	
Carbon Content	≥ 95 %
Width	Confirmed
100 mm	
Length	Confirmed
100 mm	
Thickness	Confirmed
1-5 um	

Specification: PRD.0.ZQ5.10000085348

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.