


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium 2-oxobutyrate - powder

Product Number: K0875
CAS Number: 2013-26-5
MDL: MFCD00064194
Formula: C₄H₅NaO₃
Formula Weight: 124.07 g/mol
Storage Temperature: 2 - 8 °C



| TEST | Specification |
|-------------------------------|-----------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder |
| Solubility (Color) | Colorless |
| Solubility (Turbidity) | Clear |
| 100 mg/mL in H ₂ O | |
| Enzymatic Purity (anhydrous) | ≥ 97 % |
| ¹ H NMR Spectrum | Conforms to Structure |

Specification: PRD.1.ZQ5.10000040488

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.