


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(±)-3-Amino-1,2-propanediol - 97%

Product Number: A76001
CAS Number: 616-30-8
MDL: MFCD00008140
Formula: C₃H₉NO₂
Formula Weight: 91.11 g/mol



TEST

Specification

Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid or Viscous Liquid
Infrared spectrum	Conforms to Structure
Color Test	≤ 300 APHA
Purity (Titration by HCL04)	96.5 - 103.5 %
Purity GC VPCB	≥ 96.5 %

Specification: PRD.0.ZQ5.10000022482

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.