

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cytidine 2':3'-cyclic monophosphate monosodium salt - ≥95% (HPLC)

Product Number: **C9630**
CAS Number: 15718-51-1
Formula: C₉H₁₁N₃NaO₇P
Formula Weight: 327.2 AMU



TEST

Specification

Appearance (Colour)	White
Appearance (Form)	Powder
Purity (HPLC)	≥ 95 %
Solubility (Colour)	Colorless
Solubility (Turbidity)	Clear
Water	≤ 12 %
ICP Major Analysis	Corresponds to Requirements
Identity by ICP - Sodium	Positive

Specification: PRE.0.ZQ5.10000021057

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.