

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Agar - plant cell culture tested, powder

Product Number:

**A7921**

CAS Number:

9002-18-0

MDL:

MFCD00081288

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Loss on Drying	≤ 10.0 %
Plant Cell Culture Test	Pass
pH	5.8 - 7.5
1.5% in H <sub>2</sub> O	
Residue on ignition (Ash)	≤ 2 %
Congealing Temperature	30 - 35 °C

Specification: PRD.0.ZQ5.10000015246

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.