


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Glycine anhydride – cyclic dipeptide

Product Number: G7251
CAS Number: 106-57-0
MDL: MFCD00006009
Formula: C₄H₆N₂O₂
Formula Weight: 114.10 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
at 20 mg/ml in H ₂ O	
Infrared spectrum	Conforms to Structure
Carbon	41.3 - 43.0 %
Nitrogen	24.1 - 25.0 %
Purity (TLC)	≥ 99 %
Recommended Retest Period:	-----
5 Years	

Specification: PRD.0.ZQ5.10000035762

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.