

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Pituitary Adenylate Cyclase Activating Polypeptide Type 1 Receptor antibody produced in rabbit - affinity isolated antibody, buffered aqueous solution

Product Number: **P8872**

MDL: MFCD06798507

Storage Temperature: -20 °C

TEST	Specification
Formulation Solution in phosphate buffered saline containing 0.1% sodium azide	Conforms
Working Dilution Immunohistochemistry: a working concentration of > or = 7 mg/mL is recommended using human pituitary tissues.	Conforms
Protein Content 1 mg/ml	Conforms

Specification: PRD.1.ZQ5.10000040776

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.