

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Pig IgG (whole molecule)-FITC antibody produced in rabbit - IgG fraction of antiserum, buffered aqueous solution

Product Number:

F4762

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	10.0 - 20.0 mg/ml
F/P Molar Ratio	2.5 - 6.5
IF Cells	Pig Spleen or Pig Peripheral Blood Lymphocytes
IF-Titer	1:32-1:64
Recommended Retest Period	-----
3 Years	

Specification: PRE.2.ZQ5.10000017172

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.