


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
11-Mercaptoundecanoic acid - 98%

Product Number: 674427
CAS Number: 71310-21-9
MDL: MFCD00022096
Formula: C₁₁H₂₂O₂S
Formula Weight: 218.36 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Conforms to Requirements
Powder, Crystals or Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %
Carbon	59.0 - 62.0 %
Sulfur (S)	14.3 - 15.0 %

Specification: PRD.0.ZQ5.10000054852

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.