

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris(2,4-di-tert-butylphenyl) phosphite - 98%

Product Number: 441791
CAS Number: 31570-04-4
MDL: MFCD00072706
Formula: C₄₂H₆₃O₃P
Formula Weight: 646.92 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder, Crystals, Crystalline Powder or Granules	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	76.0 - 79.9 %
Solubility (Turbidity) c = 25 mg/mL; CHCl ₃	Clear
Solubility (Color)	Colorless

Specification: PRD.1.ZQ5.10000025624

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.