

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

ProteoMass™ Apomyoglobin MALDI-MS Standard - vial of 10 nmol, (M+H⁺) 16,952.27 Da by calculation

Product Number:

A8971

CAS Number:

9008-45-1

MDL:

MFCD00198077

TEST	Specification
------	---------------

Protein (nanomoles/vial)

8.0 - 18.0

µg/vial

135.6 - 305.1

Suitability

Suitable

Suitable for use as a standard for MALDI Mass Spectrometry.

Specification Date : 04/05/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.