

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Magnesium sulfate – anhydrous, reagent grade, ≥97%

Product Number:

208094

CAS Number:

7487-88-9

MDL:

MFCD00011110

Formula:

MgSO<sub>4</sub>

Formula Weight:

120.37 g/mol

MgSO<sub>4</sub>**TEST****Specification**

Appearance (Color)

Conforms to Requirements

White to Off-White (May contain grey specks)

Appearance (Form)

Granules

X-Ray Diffraction

Conforms to Structure

Compleximetric EDTA

19.6 - 20.8 %

% Magnesium

Specification: PRD.0.ZQ5.10000036646

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.