


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(±)-3-Chloro-1,2-propanediol - 98%

Product Number: 107271
CAS Number: 96-24-2
MDL: MFCD00004712
Formula: C₃H₇ClO₂
Formula Weight: 110.54 g/mol



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃ After Oxygen Combustion	97.5 - 102.5 %
Purity (GC)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000026818

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.