


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
6-Maleimidohexanoic acid N-hydroxysuccinimide ester - ≥98%, powder

Product Number: M9794
CAS Number: 55750-63-5
MDL: MFCD00043043
Formula: C₁₄H₁₆N₂O₆
Formula Weight: 308.29 g/mol
Storage Temperature: -20 °C



TEST	Specification
Solubility (Color)	Colorless
Solubility (Turbidity) 10%, CHCl ₃	Clear
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
% Purity (HPLC)	≥ 98
Carbon	53.2 - 55.9 %
Nitrogen	8.8 - 9.4 %

Specification: PRD.0.ZQ5.10000049867

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.