


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Amino-2,6-dichloropyridine - 97%

Product Number: 565342
CAS Number: 2587-02-2
MDL: MFCD00052832
Formula: C₅H₄Cl₂N₂
Formula Weight: 163.00 g/mol



TEST

Specification

Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	96.5 - 103.5 %
After Oxygen Combustion	
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000014804

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.