

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
Quercetin - $\geq 95\%$ (HPLC), solid

Product Number: **Q4951**

CAS Number: 117-39-5

Formula: C₁₅H₁₀O₇

Formula Weight: 302.24 g/mol



TEST	Specification
Appearance (Color) Yellow	Conforms
Appearance (Form)	Powder
¹ H NMR Spectrum	Conforms to Structure
Loss on Drying	$\leq 4\%$
Purity (HPLC)	$\geq 95\%$

Specification: PRD.4.ZQ5.10000025758

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.