

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Cyclodextrin, sulfated sodium salt - extent of labeling: 12-15 mol per mol β -CD

Product Number:


389153

CAS Number:

37191-69-8

MDL:

MFCD00197857



| TEST | Specification |
|--|-----------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder or Chunks |
| Infrared spectrum | Conforms to Structure |
| Optical Rotation | |
| Typically + 45 to + 85 Degrees, c = 1%, H ₂ O | |
| Loss on Drying | |
| Typically < 20%, 2 Hours, 82 Degrees Celsius | |

Specification: PRD.1.ZQ5.10000009931

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.